



B A: Body and rgan Analyzer

- [BOA: Body and Organ Analyzer](#)
 - [What is it?](#)
 - [How does it work?](#)
 - [How to run?](#)
 - [Environment Variables](#)
 - [Notes on RabbitMQ](#)
 - [Run](#)
 - [Send a study to the analyzer](#)
 - [Notes on Performance](#)
 - [Outputs](#)
 - [Command Line Tool](#)

What is it?

BOA is a tool for segmentation of CT scans developed by the [SHIP-AI group at the Institute for Artificial Intelligence in Medicine](#). Combining the [TotalSegmentator](#) and the [Body Composition Analysis](#), this tool is capable of analyzing medical images and identifying the different structures within the human body, including bones, muscles, organs, and blood vessels.

If you use this tool, please make sure to cite the following papers: [BCA](#), [TotalSegmentator](#) and [nnU-Net](#).

How does it work?

- The user sends a study to the analyzer.
- The study is received by the Orthanc instance, which then creates a task.
- The task is picked up by a task management system, which then starts the segmentations and the computation of the Excel file with the measurements.
- If specified, the segmentations are saved locally and uploaded to the DicomWeb instance.

- If desired, we can provide an additional system where the segmentations can be viewed. Just write us an email.

are not specified, the segmentations will be lost. **Note:** Please create the directory before if you are not using the root user, this avoids having to change the permissions later.

DICOMWEB_UPLOAD_DIRECTORY : Full link to a DicomWeb instance to upload the segmentations, e.g. `https://192.168.1.4/orthanc/segmentations/`. If not specified, the segmentations are only saved locally. In this case, it is important to not use our already running Orthanc instance as the target for storing the segmentations. However, another Orthanc instance can be used.

- AccessionNumber_SeriesNumber_SeriesDescription.xlsx
- total.nii.gz
- body-regions.nii.gz
- ...

The `FP_IN_TN` option has the drawback that:

- The folder structure is not unique if your DICOMs are anonymized, i.e. if you have two patients with the same anonymized name and birthdate (e.g. John Doe, 1970), then you will have clashes.
- The output folders will not be anonymized and will contain the patient name and birthdate.

Notes on RabbitMQ

RabbitMQ always checks whether the tasks have been received by the consumer, and if the consumer takes too long, the tasks are going to be killed. This is a problem in our case, because the tasks may take more time to complete, and if 40-50 studies have been sent at the same time, they are going to be killed by RabbitMQ because they have not yet been completed. This can be set with the `—` variable of RabbitMQ, and since this highly depends on the amount of studies that are being sent in one go, we decided to disable this variable. This is done by creating a file in `/ / /` with the following content.

```
%%
[
    [ ,
    ] -
]
```

If you are using your own RabbitMQ instance, please also use this setting, or contact us if you have a better idea!

Run

Load the docker images:

```
P # !
```

or clone the repository and build the images

```
—
```

Download the [docker-compose.yml](#) file and run the following command

```

- -

```

with `python3 main.py --local-gpu` if you want to use a local GPU and `python3 main.py --triton` if you have Triton instance running. Remove `--triton` in case you already have an instance running.

You can also just

```

-

```

!!!!IMPORTANT!!!!: if you are using windows, substitute the `python3` with `python`. `python3` has not been tested extensively so if you have any problems please contact us!

Send a study to the analyzer

You can then add the instance to your PACS of choice by adding `127.0.0.1` and the port `4444` to the location manager to your PACS. Below there is a screenshot of how this looks in Horos.

Example in Horos:

⊙	Address	AETitle	Port	Q...	Retrieve	Send	TLS	Name	Send Transfer Syntax
✓	127.0.0.1	Horos	4444	<input type="checkbox"/>	C-MOVE ⇅	✓	No	⊕ This is an example	Explicit Little Endian ⇅
✓	127.0.0.1	BOA	4242	✓	C-MOVE ⇅	✓	No	⊕ Orthanc for BOA	Explicit Little Endian ⇅
✓	127.0.0.1	Onco	4242	✓	C-MOVE ⇅	<input type="checkbox"/>	No	⊕ Oncology Endpoint	Explicit Little Endian ⇅
✓	127.0.0.1	SegViewer	4243	✓	C-MOVE ⇅	✓	No	⊕ Viewer	Explicit Little Endian ⇅

In this case, the IP is the same as the one of my machine because I am testing locally. The AETTitle that you specify will be the name of the folder where the results will be stored, so you that can create different endpoints to computed different cohorts.

Notes on Performance

To make an estimate on how much power and time is needed to process a study, we used the following table provided by the TotalSegmentator. However, for very large series (e.g. 1600 slices 1mm), the performance may be worse and more CPU power may be needed. According to our tests, 16GB of GPU should be absolutely sufficient.

Outputs

The outputs of BOA are listed below, all of them will appear in the `BOA_NM_AET` folder, some of them will be uploaded to SMB and some will be uploaded with DicomWeb, if they are configured. Currently, the produced DICOM-segs have placeholders for the anatomical names of the tissues.

- Segmentations (BCA and TotalSegmentator), present in `BOA_NM_AET` and uploaded to `0` DicomWeb (optional).
 - `Total Body Segmentation ()`: Segmentation of 104 body regions (`TotalSegmentator`).
 - Intracerebral Hemorrhage Segmentation ().

- Lung Vessels and Airways Segmentation (): Segmentation of trachea/bronchia/airways (paper).
- Liver Vessels and Tumor Segmentation ().
- Hip Implant Segmentation ().
- Coronary Arteries Segmentation (): Segmentation of the coronary arteries.
- Pleural Pericardial Effusion Segmentation (): pleural effusion (paper), pericardial effusion (cite paper).
- Body Regions Segmentation (): Segmentation of the body regions (BCA).
- Body Parts Segmentation (): Segmentation of body and extremities.
- Tissue Segmentation (): Segmentation of the tissues (BCA).
- Measurements/Reports, present in and uploaded to SMB (optional).
 - : Excel file with all the measurements from the BCA and the TotalSegmentator. The sheet contains general information about the patient, such as IDs, dates, DICOM tags, contrast phase, the BOA version. The sheet has information about the segmentations that were computed together with their volume and some other statistical information. contains all the measurements for the aggregated regions of the BCA, which are also visible in the report. contains information about the volume of each tissue for each slice of the CT scan. contains the same information, but the extremities are removed from the computation.
 - : Report of the BCA findings.
 - : Preview of the TotalSegmentator segmentation.
- Other files containing measurements are also stored in the output directory (in format), and are not uploaded anywhere else. These measurements all appear in the resulting Excel report.

Command Line Tool

Additionally, the BOA can also be run from the command line to compute all the segmentations in one go and without connecting to the PACS.

First, get the image.

```
P # !
```

or clone the repository and build the image

```
- - / . - -- . /
```

then you can run your image!

Run on Linux

-- \
 .: I E L F N / T
 I: -I W G_ N / \ O
 LH ELI :P _W G A _ / T A OT \ # A

 -- \ =
 -- \
 : --\$ - - D #) IT)
 -- - -- 16 708-6 - -- 1 = 4 \ =
 -- / / \
 - / - \
 - \

 / / -- " -- . . -- - / -- -
 / / -- " --

[illegible]

- In `--models` you can specify which models you want to run: `--models seg_01` for the TotalSegmentator, `--models seg_02` for the Body Composition Analysis, `--models seg_03` for both or `--models seg_04` for all the possible models: `--models seg_01,seg_02,seg_03,seg_04`.
- You can also specify whether you want to extract the radiomics features by adding `--radiomics`.
CAREFUL: This has currently not been tested extensively, so it might not work as expected.
- During the process some segmentations are generated, which are then postprocessed and the original segmentations are deleted. If you want all versions of the segmentations, you can add `--keep-segmentations`.

- There are other parameters that either belong to the BCA or to the TotalSegmentator, which you can view in [this link](#).